

Name:
Group: GM, DM

School year:
Date of measurement:

Report

Topic: Audiometry

Exercise 1: Measure the hearing threshold of both of your ears, using air conduction. Evaluate these measurements graphically and construct the audiograms.

Devices and implements: PC-Audiometer Maico MA33

Procedure: according to the instructions.

Measured values :

Table 1:

frequency [Hz]		125	250	500	1000	2000	4000	8000
intensity level [dB]	left ear							
	right ear							

Exercise 2: Masking. Measure the masking threshold for no tested ear.

Measured values: Set the masking intensity to 25 dB

Table 2:

frequency [Hz]		125	250	500	1000	2000	4000	8000
intensity level [dB]	left ear							
	right ear							

Conclusions and commentary: